

Women's Health: Osteoporosis Screening Rule

Rule: Conceptual Structure

Contract: VA118-16-D-1008, Task Order (TO): VA-118-16-F-1008-0007, CLIN0007BA

Department of Veterans Affairs (VA)



**Knowledge Based Systems (KBS)
Office of Informatics and Information Governance (OIIG)
Clinical Decision Support (CDS)**

Publication date 06/05/2018

Version: 1.0

Women's Health: Osteoporosis Screening Rule: Rule: Conceptual Structure

by Knowledge Based Systems (KBS), Office of Informatics and Information Governance (OIIG), and Clinical Decision Support (CDS)

Publication date 06/05/2018

Copyright © 2018 B3 Group, Inc.

Copyright © 2018 Cognitive Medical Systems, Inc.

B3 Group, Inc.

NOTICE OF GOVERNMENT COPYRIGHT LICENSE AND UNLIMITED RIGHTS LICENSE

Licensed under the Apache License, Version 2.0 (the "License"); you may not use this file except in compliance with the License.

You may obtain a copy of the License at <http://www.apache.org/licenses/LICENSE-2.0>

Unless required by applicable law or agreed to in writing, software distributed under the License is distributed on an "AS IS" BASIS, WITHOUT WARRANTIES OR CONDITIONS OF ANY KIND, either express or implied. See the License for the specific language governing permissions and limitations under the License.

Portions of this content are derivative works from content produced by Cognitive Medical Systems, Inc. licensed under the Apache License, Version 2.0.

Additional portions of this content are derivative works from content contributed by Motive Medical Intelligence Inc., under Creative Commons Attribution-ShareAlike 4.0.

Contributions from 2013-2018 were performed either by US Government employees, or under US Veterans Health Administration contracts.

US Veterans Health Administration contributions by government employees are work of the U.S. Government and are not subject to copyright protection in the United States. Portions contributed by government employees are USGovWork (17USC §105). Not subject to copyright.

See: <https://www.usa.gov/government-works>

Contribution by contractors to the US Veterans Health Administration during this period are contractually contributed under the Apache License, Version 2.0 and US Government sponsorship is acknowledged under Contract VA118-16-D-1008, Task Order VA11817F10080007.

Cognitive Medical Systems, Inc.

Licensed under the Apache License, Version 2.0 (the "License"); you may not use this file except in compliance with the License.

You may obtain a copy of the License at <http://www.apache.org/licenses/LICENSE-2.0>

Unless required by applicable law or agreed to in writing, software distributed under the License is distributed on an "AS IS" BASIS, WITHOUT WARRANTIES OR CONDITIONS OF ANY KIND, either express or implied. See the License for the specific language governing permissions and limitations under the License.

This and related content produced by Cognitive Medical Systems, Inc. licensed under the Apache License, Version 2.0 is available at: <https://bitbucket.org/cogmedsys/hl7-kas-examples>

Additional portions of this content are derivative works from content contributed by Motive Medical Intelligence Inc., under Creative Commons Attribution-ShareAlike 4.0. <https://bitbucket.org/cogmedsys/kas-source-material>

Contributions from 2013-2018 were performed either by US Government employees, or under US Veterans Health Administration contracts.

US Veterans Health Administration contributions by government employees are work of the U.S. Government and are not subject to copyright protection in the United States. Portions contributed by government employees are USGovWork (17USC §105). Not subject to copyright. See: <https://www.usa.gov/government-works>

Contribution by contractors to the US Veterans Health Administration during this period are contractually contributed under the Apache License, Version 2.0 and US Government sponsorship is acknowledged under Contract VA118-16-D-1008-0007.

Table of Contents

Preface	v
Artifact Applicability	vi
Models	vii
1. External Data Definitions	1
Definitions	1
Triggers	3
2. Expression Definitions	5
3. Actions	6
4. Tabular List	10
5. Behavior Symbols	13
A. References	15

List of Tables

1. Revision History	v
2. Clinical White Paper Contributors	v
3. Artifact Identifier	v
4. Applicability Foci, Description and Codes	vi
5. Model References	vii
1.1. PatientAge	1
1.2. DiagnosedWithOsteoporosis	1
1.3. DXAScanning	1
1.4. OsteoporosisScreeningInterval	1
1.5. Postmenopausal	1
1.6. HasHistoryOfFragilityFracture	1
1.7. LowBodyWeight	2
1.8. BoneLossCausedByMedication	2
1.9. BoneLossCausedByDisease	2
1.10. ParentalMedicalHistoryOfHipFracture	2
1.11. CurrentTobaccoSmoker	2
1.12. HasAlcoholism	2
1.13. RheumatoidArthritis	2
1.14. DXAScanTScore	3
1.15. MostRecentDXAScanTScore	3
1.16. AccessOfPatientRecord	3
1.17. DxaScanGreaterEqualFifteenYears	3
1.18. DxaScanGreaterEqualFiveYears	3
1.19. DxaScanGreaterEqualOneYear	3
1.20. EndPatientPerferredScreeningInterval	4
4.1. Terminology Versions	10
4.2. Terminology References	10
5.1. Group Organizational Behavior	13
5.2. Group Selection Behavior	13
5.3. Required Behavior	13
5.4. Precheck Behavior	14
5.5. Cardinality Behavior	14
5.6. Item Flags	14
5.7. Read Only Behavior	14

Preface

Table 1. Revision History

Date	Life Cycle Event
June 5, 2018	Published
June 5, 2018	Reviewed
April 25, 2018	Published
April 25, 2018	Reviewed
March 16, 2018	Created

Table 2. Clinical White Paper Contributors

Name	Role	Affiliation
C. Yvette Williams-Harris, MD MPH	Author	Physician, 1670 Clairmont Rd, Atlanta, GA 30033
Megan Gerber, MD MPH	Author	Physician, VA Boston Health Care System, 150 South Huntington Ave, Boston, MA 02130
Caitlin M. Cusack, MD MPH	Author	Contractor-Physician Informaticist, Office of Women's Health Services, 8455 Colesville Rd, 12th Floor, Silver Spring, MD 20910

Table 3. Artifact Identifier

Domain	Artifact ID	Name
urn:va.gov:kbs:knart:artifact:rl	6ab2870a-7bc5-53c9-be4e-9fb798c50e15	O8

Artifact Applicability

Table 4. Applicability Foci, Description and Codes

Focus	Description	Code System	Code	Value Set	Value Set Version
PatientGender	Female	SNOMED CT	248152002 Female (finding)	N/A	N/A
PatientAgeGroup	Population 18 years old or older	SNOMED CT	133936004 Adult (person)	N/A	N/A
WorkflowSetting	Routine			N/A	N/A
TargetUser	Primary Care	SNOMED CT	453231000124104 Primary care provider (occupation)	N/A	N/A
TargetUser	Gynecology	SNOMED CT	702892004 Obstetrics and gynecology clinic (environment)	N/A	N/A
ClinicalVenue	Outpatient	SNOMED CT	33022008 Hospital-based outpatient department (environment)	N/A	N/A

Models

Table 5. Model References

Referenced Model	Description
urn:solor.io:anf-model:1.0	VA Analysis Normal Form Model

Chapter 1. External Data Definitions

Definitions

Table 1.1. PatientAge

Expression: type=elm:SingletonFrom Annotation: Codes: elm:value[elm:Code]: [398166005 Performed (qualifier value)]elm:value[elm:Code]: [[a997cc03-3e99-40eb-833a-6374c7750a3a Observation procedure (procedure)] ->(363702006 Has focus (attribute))->[105726004 Age AND/OR growth finding (finding)]]
--

Table 1.2. DiagnosedWithOsteoporosis

Expression: type=elm:Query Annotation: Codes: elm:value[elm:Code]: [398166005 Performed (qualifier value)]elm:value[elm:Code]: [[a997cc03-3e99-40eb-833a-6374c7750a3a Observation procedure (procedure)] ->(363702006 Has focus (attribute))->[64859006 Osteoporosis (disorder)]]

Table 1.3. DXAScanning

Expression: type=elm:Query Annotation: Codes: elm:value[elm:Code]: [398166005 Performed (qualifier value)]elm:value[elm:Code]: [241686001 Dual energy X-ray photon absorptiometry (procedure)]

Table 1.4. OsteoporosisScreeningInterval

Expression: type=elm:Query Annotation: Codes: elm:value[elm:Code]: [398166005 Performed (qualifier value)]elm:value[elm:Code]: [300004007 Screening for osteoporosis (procedure)]
--

Table 1.5. Postmenopausal

Expression: type=elm:Query Annotation: Codes: elm:value[elm:Code]: [398166005 Performed (qualifier value)]elm:value[elm:Code]: [[a997cc03-3e99-40eb-833a-6374c7750a3a Observation procedure (procedure)] ->(363702006 Has focus (attribute))->[76498008 Postmenopausal state (finding)]]
--

Table 1.6. HasHistoryOfFragilityFracture

Expression: type=elm:Query Annotation: Codes: elm:value[elm:Code]: [398166005 Performed (qualifier value)]elm:value[elm:Code]: [[a997cc03-3e99-40eb-833a-6374c7750a3a Observation procedure (procedure)] ->(363702006 Has focus (attribute))->[268029009 Pathological fracture (disorder)]]

Table 1.7. LowBodyWeight

Expression: type=elm:Query
Annotation:
Codes: elm:value[elm:Code]: [398166005 Performed (qualifier value)]elm:value[elm:Code]: [248342006 Underweight (finding)]

Table 1.8. BoneLossCausedByMedication

Expression: type=elm:Query
Annotation:
Codes: elm:value[elm:Code]: [398166005 Performed (qualifier value)]elm:value[elm:Code]: [[a997cc03-3e99-40eb-833a-6374c7750a3a Observation procedure (procedure)] ->(363702006 Has focus (attribute))>[14651005 Drug-induced osteoporosis (disorder)]]

Table 1.9. BoneLossCausedByDisease

Expression: type=elm:Query
Annotation:
Codes: elm:value[elm:Code]: [398166005 Performed (qualifier value)]elm:value[elm:Code]: [[a997cc03-3e99-40eb-833a-6374c7750a3a Observation procedure (procedure)] ->(363702006 Has focus (attribute))>[17b00578-6b1b-4e0b-b2ee-76f69e6061ec Osteoporosis due to disease (disorder)]]

Table 1.10. ParentalMedicalHistoryOfHipFracture

Expression: type=elm:Query
Annotation:
Codes: elm:value[elm:Code]: [398166005 Performed (qualifier value)]elm:value[elm:Code]: [[a997cc03-3e99-40eb-833a-6374c7750a3a Observation procedure (procedure)] ->(363702006 Has focus (attribute))>[263225007 Fracture of proximal end of femur (disorder)]]

Table 1.11. CurrentTobaccoSmoker

Expression: type=elm:Query
Annotation:
Codes: elm:value[elm:Code]: [398166005 Performed (qualifier value)]elm:value[elm:Code]: [[a997cc03-3e99-40eb-833a-6374c7750a3a Observation procedure (procedure)] ->(363702006 Has focus (attribute))>[77176002 Smoker (finding)]]

Table 1.12. HasAlcoholism

Expression: type=elm:Query
Annotation:
Codes: elm:value[elm:Code]: [398166005 Performed (qualifier value)]elm:value[elm:Code]: [[a997cc03-3e99-40eb-833a-6374c7750a3a Observation procedure (procedure)] ->(363702006 Has focus (attribute))>[7200002 Alcoholism (disorder)]]

Table 1.13. RheumatoidArthritis

Expression: type=elm:Query

Annotation:

Codes: elm:value[elm:Code]: [398166005 |Performed (qualifier value)|elm:value[elm:Code]:
[[a997cc03-3e99-40eb-833a-6374c7750a3a |Observation procedure (procedure)] ->(363702006 |Has
focus (attribute))>[69896004 |Rheumatoid arthritis (disorder)]]

Table 1.14. DXAScanTScore

Expression: type=elm:Query

Annotation:

Codes: elm:value[elm:Code]: [398166005 |Performed (qualifier value)|elm:value[elm:Code]:
[[a997cc03-3e99-40eb-833a-6374c7750a3a |Observation procedure (procedure)] ->(363702006 |Has
focus (attribute))>[3dfd49ed-a7e4-41a5-ab44-356015e345a2 |Dual energy X-ray photon absorptiome-
try scan T score]]

Table 1.15. MostRecentDXAScanTScore

Expression: type=elm:Indexer

Annotation:

Codes:

Triggers

Table 1.16. AccessOfPatientRecord

Trigger: type=DataEventTrigger, DataElementAccessed

Expression: type=elm:Instance

Annotation: Access of the patient record in a preventive care encounter at a primary care clinic or in an
Obstetrics/Gynecology (OB/GYN) or women's health outpatient department.

Codes: elm:value[elm:Code]: [398166005 |Performed (qualifier value)|elm:value[elm:Code]: [TSR-
NoCode]

Table 1.17. DxaScanGreaterEqualFifteenYears

Trigger: type=DataEventTrigger, DataElementAccessed

Expression: type=elm:Instance

Annotation: Previous DXA scan >= 15 years before today.

Codes: elm:value[elm:Code]: [398166005 |Performed (qualifier value)|elm:value[elm:Code]:
[241686001 |Dual energy X-ray photon absorptiometry (procedure)]]

Table 1.18. DxaScanGreaterEqualFiveYears

Trigger: type=DataEventTrigger, DataElementAdded

Expression: type=elm:Instance

Annotation: Previous DXA scan >= 5 years before today.

Codes: elm:value[elm:Code]: [398166005 |Performed (qualifier value)|elm:value[elm:Code]:
[241686001 |Dual energy X-ray photon absorptiometry (procedure)]]

Table 1.19. DxaScanGreaterEqualOneYear

Trigger: type=DataEventTrigger, DataElementAdded

Expression: type=elm:Instance
 Annotation: Previous DXA scan >= 1 year before today.
 Codes: elm:value[elm:Code]: [398166005 |Performed (qualifier value)]elm:value[elm:Code]:
 [241686001 |Dual energy X-ray photon absorptiometry (procedure)]

Table 1.20. EndPatientPerferredScreeningInterval

Trigger: type=DataEventTrigger, DataElementModified
 Expression: type=elm:Instance
 Annotation: End of set patient-preferred screening interval.
 Codes: elm:value[elm:Code]: [385644000 |Requested (qualifier value)]elm:value[elm:Code]:
 [300004007 |Screening for osteoporosis (procedure)]

Chapter 2. Expression Definitions

No expression definitions are present.

Chapter 3. Actions

Condition: elm:GreaterOrEqual (PatientAge elm:Quantity(65)) Condition: elm:Not (elm:Exists(DiagnosedWithOsteoporosis)) Condition: elm:Not (elm:Exists(DXAScanning)) Condition: elm:Not (elm:Exists(OsteoporosisScreeningInterval))
Outpatient Visit, >= 65 Years of Age, No Prior Dual-Energy X-ray Absorptiometry (DXA)
☆ Identify the patient as an osteoporosis screening candidate. actionSentence[type=elm:Instance, classType=anf:ClinicalStatement] <div> statementType: Precoordinated Expression 385644000 Re- requested (qualifier value) topic: Precoordinated Expression TSR-NoCode </div> ☆ Open Documentation Template: Osteoporosis Screening KNART. actionSentence[type=elm:Instance, classType=anf:ClinicalStatement] <div> statementType: Precoordinated Expression 385644000 Re- requested (qualifier value) topic: Precoordinated Expression TSR-NoCode </div> ☆ Open Order Set: Osteoporosis Screening KNART. actionSentence[type=elm:Instance, classType=anf:ClinicalStatement] <div> statementType: Precoordinated Expression 385644000 Re- requested (qualifier value) topic: Precoordinated Expression TSR-NoCode </div> ☆ Make screening guideline recommendations available for review. actionSentence[type=elm:Instance, classType=anf:ClinicalStatement] <div> statementType: Precoordinated Expression 385644000 Re- requested (qualifier value) topic: Precoordinated Expression TSR-NoCode </div>
Condition: elm:Less (PatientAge elm:Quantity(65)) Condition: elm:Exists (Postmenopausal) Condition: elm:AnyTrue Condition: elm:Not (elm:Exists(DiagnosedWithOsteoporosis)) Condition: elm:Not (elm:Exists(OsteoporosisScreeningInterval))
Outpatient Visit, < 65 Years of Age
☆ Identify the patient as an osteoporosis screening candidate actionSentence[type=elm:Instance, classType=anf:ClinicalStatement] <div> statementType: Precoordinated Expression 385644000 Re- requested (qualifier value) topic: Precoordinated Expression TSR-NoCode </div> ☆ Open Documentation Template: Osteoporosis Screening KNART actionSentence[type=elm:Instance, classType=anf:ClinicalStatement]

statementType: Precoordinated Expression 385644000 |Requested (qualifier value)|
 topic: Precoordinated Expression TSR-NoCode

☆ Open Order Set: Osteoporosis Screening KNART

actionSentence[type=elm:Instance, classType=anf:ClinicalStatement]

statementType: Precoordinated Expression 385644000 |Requested (qualifier value)|
 topic: Precoordinated Expression TSR-NoCode

☆ Make screening guideline recommendations available for review

actionSentence[type=elm:Instance, classType=anf:ClinicalStatement]

statementType: Precoordinated Expression 385644000 |Requested (qualifier value)|
 topic: Precoordinated Expression TSR-NoCode

Condition:elm:GreaterOrEqual (PatientAge elm:Quantity(65))
Condition:elm:GreaterOrEqual (MostRecentDXAScanTScore elm:Quantity(-1.5))
Condition:elm:Not (elm:Exists(DiagnosedWithOsteoporosis))
Condition:elm:Not (elm:Exists(OsteoporosisScreeningInterval))

<= 15 Years Since Previous Dual-Energy X-ray Absorptiometry (DXA) Scan, >= 65 Years of Age, Most Recent T-Score >= -1.5

☆ Identify the patient as an osteoporosis screening candidate

actionSentence[type=elm:Instance, classType=anf:ClinicalStatement]

statementType: Precoordinated Expression 385644000 |Requested (qualifier value)|
 topic: Precoordinated Expression TSR-NoCode

☆ Open Documentation Template: Osteoporosis Screening KNART

actionSentence[type=elm:Instance, classType=anf:ClinicalStatement]

statementType: Precoordinated Expression 385644000 |Requested (qualifier value)|
 topic: Precoordinated Expression TSR-NoCode

☆ Open Order Set: Osteoporosis Screening KNART

actionSentence[type=elm:Instance, classType=anf:ClinicalStatement]

statementType: Precoordinated Expression 385644000 |Requested (qualifier value)|
 topic: Precoordinated Expression TSR-NoCode

☆ Make screening guideline recommendations available for review

actionSentence[type=elm:Instance, classType=anf:ClinicalStatement]

statementType: Precoordinated Expression 385644000 |Requested (qualifier value)|
 topic: Precoordinated Expression TSR-NoCode

Condition:elm:GreaterOrEqual (PatientAge elm:Quantity(65))
Condition:elm:And (elm:Greater(MostRecentDXAScanTScore elm:Quantity(-2.0))
 elm:Less(MostRecentDXAScanTScore elm:Quantity(-1.5)))

Condition: elm:Not (elm:Exists(DiagnosedWithOsteoporosis)) Condition: elm:Not (elm:Exists(OsteoporosisScreeningInterval))
5 Years Since Previous Dual-Energy X-ray Absorptiometry (DXA) Scan, >= 65 Years of Age, Most Recent T-Score is > -2.0 and < -1.5
☆ Identify the patient as an osteoporosis screening candidate actionSentence[type=elm:Instance, classType=anf:ClinicalStatement] <div>statementType: Precoordinated Expression 385644000 Re- requested (qualifier value) </div> <div>topic: Precoordinated Expression TSR-NoCode</div> ☆ Open Documentation Template: Osteoporosis Screening KNART actionSentence[type=elm:Instance, classType=anf:ClinicalStatement] <div>statementType: Precoordinated Expression 385644000 Re- requested (qualifier value) </div> <div>topic: Precoordinated Expression TSR-NoCode</div> ☆ Open Order Set: Osteoporosis Screening KNART actionSentence[type=elm:Instance, classType=anf:ClinicalStatement] <div>statementType: Precoordinated Expression 385644000 Re- requested (qualifier value) </div> <div>topic: Precoordinated Expression TSR-NoCode</div> ☆ Make screening guideline recommendations available for review actionSentence[type=elm:Instance, classType=anf:ClinicalStatement] <div>statementType: Precoordinated Expression 385644000 Re- requested (qualifier value) </div> <div>topic: Precoordinated Expression TSR-NoCode</div>
Condition: elm:GreaterOrEqual (PatientAge elm:Quantity(65)) Condition: elm:In (MostRecentDXAScanTScore elm:Interval()) Condition: elm:Not (elm:Exists(DiagnosedWithOsteoporosis)) Condition: elm:Not (elm:Exists(OsteoporosisScreeningInterval))
1 Year Since Previous Dual-Energy X-ray Absorptiometry (DXA) Scan, >= 65 Years of Age, Most Recent T-Score is >= -2.49 and <= -2
☆ Identify the patient as an osteoporosis screening candidate actionSentence[type=elm:Instance, classType=anf:ClinicalStatement] <div>statementType: Precoordinated Expression 385644000 Re- requested (qualifier value) </div> <div>topic: Precoordinated Expression TSR-NoCode</div> ☆ Open Documentation Template: Osteoporosis Screening KNART actionSentence[type=elm:Instance, classType=anf:ClinicalStatement] <div>statementType: Precoordinated Expression 385644000 Re- requested (qualifier value) </div> <div>topic: Precoordinated Expression TSR-NoCode</div> ☆ Open Order Set: Osteoporosis Screening KNART actionSentence[type=elm:Instance, classType=anf:ClinicalStatement]

	statementType: Precoordinated Expression 385644000 Requested (qualifier value)
	topic: Precoordinated Expression TSR-NoCode
☆ Make screening guideline recommendations available for review	
actionSentence[type=elm:Instance, classType=anf:ClinicalStatement]	
	statementType: Precoordinated Expression 385644000 Requested (qualifier value)
	topic: Precoordinated Expression TSR-NoCode
Condition: elm:Not (elm:Exists(DiagnosedWithOsteoporosis))	
Condition: elm:Exists (OsteoporosisScreeningInterval)	
Patient-Preferred Screening Interval	
☆ Identify the patient as an osteoporosis screening candidate	
actionSentence[type=elm:Instance, classType=anf:ClinicalStatement]	
	statementType: Precoordinated Expression 385644000 Requested (qualifier value)
	topic: Precoordinated Expression TSR-NoCode
☆ Open Documentation Template: Osteoporosis Screening KNART	
actionSentence[type=elm:Instance, classType=anf:ClinicalStatement]	
	statementType: Precoordinated Expression 385644000 Requested (qualifier value)
	topic: Precoordinated Expression TSR-NoCode
☆ Open Order Set: Osteoporosis Screening KNART	
actionSentence[type=elm:Instance, classType=anf:ClinicalStatement]	
	statementType: Precoordinated Expression 385644000 Requested (qualifier value)
	topic: Precoordinated Expression TSR-NoCode
☆ Make screening guideline recommendations available for review	
actionSentence[type=elm:Instance, classType=anf:ClinicalStatement]	
	statementType: Precoordinated Expression 385644000 Requested (qualifier value)
	topic: Precoordinated Expression TSR-NoCode

Chapter 4. Tabular List

Terminology Service Request (TSR) Mappings

Table 4.1. Terminology Versions

Name	Identifier	Version
SNOMED CT	2.16.840.1.113883.6.96	United States Edition 20180301
RXNorm	2.16.840.1.113883.6.88	5-Mar-18

Table 4.2. Terminology References

System	Code	Display Text ^a	References ^b
SNOMED CT	133936004 Adult (person)	Population 18 years old or older	1
SNOMED CT	241686001 Dual energy X-ray photon absorptiometry (procedure)	Precoordinated Expression	4
SNOMED CT	248152002 Female (finding)	Female	1
SNOMED CT	248342006 Underweight (finding)	Precoordinated Expression	1
SNOMED CT	300004007 Screening for osteoporosis (procedure)	Precoordinated Expression	2
SNOMED CT	33022008 Hospital-based outpatient department (environment)	Outpatient	1
SNOMED CT	385644000 Requested (qualifier value)	Precoordinated Expression	25
SNOMED CT	398166005 Performed (qualifier value)	Precoordinated Expression	18
SNOMED CT	453231000124104 Primary care provider (occupation)	Primary Care	1
SNOMED CT	702892004 Obstetrics and gynecology clinic (environment)	Gynecology	1
SNOMED CT	TSR-NoCode ^c	Precoordinated Expression	25
SNOMED CT	[a997cc03-3e99-40e-b-833a-6374c7750a3a Observation procedure (procedure)] ->(363702006 Has focus (attribute))->[105726004 Age AND/OR growth finding (finding)]	Postcoordinated Expression	1

System	Code	Display Text ^a	References ^b
SNOMED CT	[a997cc03-3e99-40e-b-833a-6374c7750a3a Observation procedure (procedure)] ->(363702006 Has focus (attribute))->[14651005 Drug-induced osteoporosis (disorder)]	Postcoordinated Expression	1
SNOMED CT	[a997cc03-3e99-40e-b-833a-6374c7750a3a Observation procedure (procedure)] ->(363702006 Has focus (attribute))->[17b00578-6b1b-4e0b-b2ee-76f69e6061ec Osteoporosis due to disease (disorder)]	Postcoordinated Expression	1
SNOMED CT	[a997cc03-3e99-40e-b-833a-6374c7750a3a Observation procedure (procedure)] ->(363702006 Has focus (attribute))->[263225007 Fracture of proximal end of femur (disorder)]	Postcoordinated Expression	1
SNOMED CT	[a997cc03-3e99-40e-b-833a-6374c7750a3a Observation procedure (procedure)] ->(363702006 Has focus (attribute))->[268029009 Pathological fracture (disorder)]	Postcoordinated Expression	1
SNOMED CT	[a997cc03-3e99-40e-b-833a-6374c7750a3a Observation procedure (procedure)] ->(363702006 Has focus (attribute))->[3dfd49ed-a7e4-41a5-ab44-356015e345a2 Dual energy X-ray photon absorptiometry scan T score]	Postcoordinated Expression	1
SNOMED CT	[a997cc03-3e99-40e-b-833a-6374c7750a3a Observation procedure (procedure)] ->(363702006 Has focus (attribute))->[64859006 Osteoporosis (disorder)]	Postcoordinated Expression	1

System	Code	Display Text ^a	References ^b
SNOMED CT	[a997cc03-3e99-40e-b-833a-6374c7750a3a Observation procedure (procedure)] - >(363702006 Has focus (attribute))->[69896004 Rheumatoid arthritis (disorder)]	Postcoordinated Expression	1
SNOMED CT	[a997cc03-3e99-40e-b-833a-6374c7750a3a Observation procedure (procedure)] - >(363702006 Has focus (attribute))->[7200002 Alcoholism (disorder)]	Postcoordinated Expression	1
SNOMED CT	[a997cc03-3e99-40e-b-833a-6374c7750a3a Observation procedure (procedure)] - >(363702006 Has focus (attribute))->[76498008 Postmenopausal state (finding)]	Postcoordinated Expression	1
SNOMED CT	[a997cc03-3e99-40e-b-833a-6374c7750a3a Observation procedure (procedure)] - >(363702006 Has focus (attribute))->[77176002 Smoker (finding)]	Postcoordinated Expression	1

^aIf a code is used multiple times in the KNART, only the display text of the first instance is shown.

^bCount of the number of times the given code system and code pair is used in the KNART.

^cTSR-NoCode is a placeholder indicating a code was requested, but was not provided.

Chapter 5. Behavior Symbols

Table 5.1. Group Organizational Behavior

Sym- bol	Name	Definition
▶	Sentence Group	A group of related alternative actions is a sentence group if the item referenced by the action is the same in all the actions, and each action simply constitutes a different variation on how to specify the details for that item. For example, two actions that could be in a SentenceGroup are "aspirin, 500 mg, 2 times per day" and "aspirin, 300 mg, 3 times per day". In both cases, aspirin is the item referenced by the action, and the two actions represent two different options for how aspirin might be ordered for the patient. Note that a SentenceGroup would almost always have an associated selection behavior of "AtMostOne", unless it's a required action, in which case, it would be "ExactlyOne".
▷	Logical Group	A group with this behavior logically groups its sub-elements, and may be shown as a visual group to the end user, but it is not required to do so.
➤	Visual Group	Any group marked with this behavior should be displayed as a visual group to the end user.

Table 5.2. Group Selection Behavior

Sym- bol	Name	Definition
□	Any	Any number of the items in the group may be chosen, from zero to all.
⦿	All	All the items in the group must be selected as a single unit.
⊙	AllOrNone	All the items in the group are meant to be chosen as a single unit: either all must be selected by the end user, or none may be selected.
○	ExactlyOne	The end user must choose one and only one of the selectable items in the group. The user may not choose none of the items in the group.
⊛	AtMostOne	The end user may choose zero or at most one of the items in the group.
⦿	OneOrMore	The end user must choose a minimum of one, and as many additional as desired.

Table 5.3. Required Behavior

Sym- bol	Name	Definition
◆	Must	An action with this behavior must be included in the actions processed by the end user; the end user may not choose not to include this action.

Sym- bol	Name	Definition
◇	Could	An action with this behavior may be included in the set of actions processed by the end user.
➤	MustUnlessDocumented	An action with this behavior must be included in the set of actions processed by the end user, unless the end user provides documentation as to why the action was not included.

Table 5.4. Precheck Behavior

Sym- bol	Name	Definition
▲	Yes	An action with this behavior is one of the most frequent actions that is, or should be, included by an end user, for the particular context in which the action occurs. The system displaying the action to the end user should consider "pre-checking" such an action as a convenience for the user.
▽	No	An action with this behavior is one of the less frequent actions included by the end user, for the particular context in which the action occurs. The system displaying the actions to the end user would typically not "pre-check" such an action.

Table 5.5. Cardinality Behavior

Sym- bol	Name	Definition
◆	Single	An action with this behavior may only be completed once.
❖	Multiple	An action with this behavior may be repeated multiple times.

Table 5.6. Item Flags

Sym- bol	Name	Definition
☞	fillIn	This item, in a list entry, allows the user to enter a fill in value that is not present in the set of presented choices.

Table 5.7. Read Only Behavior

Sym- bol	Name	Definition
☆	true	For a particular action or action group, specifies whether the elements are read only.

Appendix A. References

This appendix contains the list of related resources and supporting documents used in creating this KNART.

List of References

Related Resources

[CCWP] *Women's Health: Osteoporosis Screening Clinical Content White Paper*

[CSD] *Women's Health: Osteoporosis Screening Rule Conceptual Structure Document*

[KVRpt] *Women's Health: Osteoporosis Screening Rule KNART Validation Report*

Supporting Evidence

[AGOG 2012] *Committee on Practice Bulletins-Gynecology, ACOG. AGOG Practice Bulletin N. 129. Osteoporosis. Obstet Gynecol. 2012;120(3):718-734 (link [<https://doi.org/10.1097/AOG.0b013e31826dc446>])*

[Cosman 2014] *Cosman F, de Beur S.J., LeBoff M.S., Lewiecki E.M., Tanner B, Randall S, Lindsay R. Clinician's Guide to Prevention and Treatment of Osteoporosis. Osteoporosis International. 2014;25(10):2359-2381 (link [<http://doi.org/10.1007/s00198-014-2794-2>])*

[Gourlay 2012] *Gourlay ML, Fine JP, Preisser JS, May RC, Li C, Lui LY, Ransohoff DF, Cauley JA, Ensrud KE. Bone-density testing interval and transition to osteoporosis in older women. New England Journal of Medicine. 2012 Jan;19(366):225-33 (link [<https://dx.doi.org/10.1056%2FNEJMoa1107142>])*

[NIH 2017] *National Institutes of Health. Calcium Fact Sheet for Health Professionals. Accessed August 2017. (link [<https://ods.od.nih.gov/factsheets/Calcium-HealthProfessional/>])*

[USPTF 2015] *USPSTF. Final Update Summary: Osteoporosis: U.S. Preventive Services Task Force website. Reviewed July 2015 (link [<https://www.uspreventiveservicestaskforce.org/Page/Document/UpdateSummaryFinal/osteoporosis-screening>])*